

Name	Use Case/Requirements	Description
power	1 and 2	User: can press this button to either turn on CES device or turn it off
ok	User is on a menu	User:The ok button allows users to traverse the menu objects on the CES device. Pressing okay while on an option will make the option load up.
up		User:The Up botton allows the user to augment the power of the 3 CES device
down		User: The Down Button will allow the user to decrease the power of the 3 CES device

timer		User: The timer button will allow the user to cycle through each interval of time option on the CES device. (20,40,and 4 60)
back	User has turned on the device and left the starting page.	User: the back button will allow the user to interact with menus as well. Pressing it returns to the previous page.
charge	battery life has been drained	Simulator: This button is used to "charge" the device, simulating either the replacement of batteries or charging

PowerSurge	<p>Simulator: This button simulates a powersurge in the device, setting it's power setting to 12 over 700.</p>
History	<p>User: The history button allows a user to access their previous sessions and 6 view them.</p>
Finish Session	<p>Session must have started</p> <p>Simulator: this button is used to simulate a whole session of CES. Using it will set the timer to zero.</p>
Change Frequency	<p>User: This button allows the User to cycle between each frequency 8 option</p>

Change Waveform		User: This button allows the User to cycle between each waveform 5 option
Realistic Power Draw		Simulator: This simulates real time battery drain. While this is active, the battery drains as it would on a normal device. If it's turned off, the battery drains much faster to test 7 the errors.
Contact	Session must be ready to start	Simulator: this check represents "contact" with the user. Meant to simulate the user using the device

User: Can  
press this  
button to  
lock their  
settings, so  
that they are  
no longer  
changed  
when other  
buttons are  
9 pushed.

Lock